

THORACIC & LUNG SURGERY



COLLIN
S-0010-19 19.0cm



STILLE
S-0020-21 21.0cm



GLUCK
S-0030-19 19.0cm
S-0030-22 22.0cm



S-0040-20 20.0cm
S-0040-23 23.0cm



S-0050-22 22.0cm



S-0060-17 17.0cm

THORACIC & LUNG SURGERY



GIERTZ -STILLE
S-0070-27 27.0cm



SAUERBRUCH - FREY
S-0080-36 36.0cm



BRUNNER
32.0cm
S-0090-28 right
S-0092-28 left



LEBSCHER
S-0100-24 24.5cm



LEBSCHER
S-0110-25 25.0cm



SAUERBRUCH
S-0120-26 26.0cm

THORACIC & LUNG SURGERY




CORYLLOS
S-0130-35 35.0cm

S-0132-35 35.0cm



BETHUNE
S-0140-34 34.0cm



DOYEN
S-0150-19 19.0cm

THORACIC & LUNG SURGERY



BETHUNE-NELSON
S-0160-34 34.0cm
with probe ended blades



ROBERTS
S-0170-34 34.0cm
with angular probe ended blades



ROBERTS
S-0180-34 34.0cm
with angular probe ended
blades of first and second ribs

THORACIC & LUNG SURGERY




SCHUMACHER
S-0190-21 21.0cm



TUDOR EDWARDS
S-0200-23 23.0cm

THORACIC & LUNG SURGERY



CASTAÑEDA				
S/S	ALU	a	b	c
S-0300-01	S-7300-01	10.5mm	32mm	65mm
S-0300-02	S-7300-02	16mm	30mm	65mm




CASTAÑEDA				
S/S	ALU	a	b	c
S-0302-01	S-7302-01	12mm	45mm	70mm



CASTAÑEDA				
S/S	ALU	a	b	c
S-0304-01	S-7304-01	15mm	60mm	120mm

THORACIC & LUNG SURGERY

Scheme of dimensions




COOLEY

S/S	ALU	a	b	c
S-0310-01	S-7310-01	12mm	13mm	95mm
S-0310-02	S-7310-02	12mm	17mm	75mm
S-0310-03	S-7310-03	14mm	16mm	85mm



COOLEY


S/S	ALU	a	b	c
S-0312-01	S-7312-01	15mm	30mm	115mm
S-0312-02	S-7312-02	15mm	30mm	140mm
S-0312-03	S-7312-03	16mm	30mm	105mm
S-0312-04	S-7312-04	20mm	45mm	140mm
S-0312-05	S-7312-05	21mm	44mm	150mm
S-0312-06	S-7312-06	30mm	35mm	140mm



COOLEY

S/S	ALU	a	b	c
S-0314-01	S-7314-01	30mm	50mm	185mm
S-0314-02	S-7314-02	35mm	50mm	200mm
S-0314-03	S-7314-03	37mm	49mm	200mm

THORACIC & LUNG SURGERY



COOLEY

for neonatal

S/S	ALU	a	b	c
S-0320-01	S-7320-01	12mm	15mm	70mm



COOLEY

infant, sternal spreader

S/S	ALU	a	b	c
S-0322-01	S-7322-01	15mm	20mm	100mm





COOLEY

for babies

S/S	ALU	a	b	c
S-0330-01	S-7330-01	20mm	25mm	100mm

THORACIC & LUNG SURGERY



COOLEY

for children

S/S	ALU	a	b	c
S-0340-01	S-7340-01	21mm	45mm	145mm

for adult

S-0340-02	S-7340-02	32mm	50mm	190mm
-----------	-----------	------	------	-------



COOLEY

for children

S/S	ALU	a	b	c
S-0350-01	S-7350-01	35mm	100mm	180mm

for adult

S-0350-02	S-7350-02	35mm	115mm	230mm
-----------	-----------	------	-------	-------



COOLEY

for adults with suture separator


S/S	ALU	a	b	c
S-0360-01	S-7360-01	35mm	100mm	180mm

THORACIC & LUNG SURGERY




H-COOLEY
with swivel convex blades, for babies

S/S	ALU	a	b	c
S-0370-01	S-7370-01	10mm	15mm	90mm




H-COOLEY
with swivel blades, for babies

S/S	ALU	a	b	c
S-0370-02	S-7370-02	10mm	15mm	90mm



STRUCK
for infant

S/S
S-0380-00
consisting of
S-0380-90 spreader only
S-0380-04 1 pair
S-0380-05 1 pair



S/S blades only

S-0380-01
10x10mm

S-0380-02
10x15mm


S-0380-03
10x20mm

S-0380-04
15x20mm

S-0380-05
20x20mm

THORACIC & LUNG SURGERY


Scheme of dimensions



WEINBERG


for children

S/S	ALU	a	b	c
S-0390-01	S-7390-01	12mm	20mm	95mm



FINOCHIETTO-INFANT


S/S	ALU	a	b	c
S-0400-01	S-7400-01	12mm	15mm	70mm
S-0400-02	S-7400-02	15mm	15mm	75mm



FINOCHIETTO-BABY-HAIGHT


S/S	ALU	a	b	c
S-0404-01	S-7404-01	27mm	30mm	95mm
S-0404-02	S-7404-02	31mm	31mm	100mm

THORACIC & LUNG SURGERY



95 mm


Scheme of dimensions



Size of blade Total spread

FINOCHIETTO-BABY

S/S	ALU	a	b	c
S-0410-01	S-7410-01	14mm	15mm	75mm
S-0410-02	S-7410-02	18mm	20mm	75mm



180 mm



S/S	ALU	a	b	c
S-0420-01	S-7420-01	28mm	32mm	120mm



FINOCHIETTO


S/S	ALU	a	b	c
S-0430-01	S-7430-01	12mm	15mm	100mm
S-0430-02	S-7430-02	28mm	30mm	143mm
S-0430-03	S-7430-03	31mm	45mm	210mm
S-0430-04	S-7430-04	42mm	61mm	245mm
S-0430-05	S-7430-05	74mm	64mm	245mm

THORACIC & LUNG SURGERY



FINOCHIETTO


S/S	ALU	a	b	c
S-0440-01	S-7440-01	36mm	45mm	140mm




FINOCHIETTO

S/S	ALU	a	b	c
S-0442-01	S-7442-01	44mm	65mm	210mm
S-0442-02	S-7442-02	65mm	65mm	210mm

THORACIC & LUNG SURGERY




Scheme of dimensions



FINOCHIETTO



S/S	ALU	a	b	c
S-0450-01	S-7450-01	30mm	45mm	140mm
S-0450-02	S-7450-02	38mm	44mm	135mm
S-0450-03	S-7450-03	52mm	62mm	175mm
S-0450-04	S-7450-04	85mm	65mm	175mm



FINOCHIETTO

S/S	ALU	a	b	c
S-0460-01	S-7460-01	40mm	40mm	165mm
S-0460-02	S-7460-02	45mm	65mm	200mm
S-0460-03	S-7460-03	75mm	65mm	250mm

THORACIC & LUNG SURGERY



FINOCHIETTO

for adults

S/S	ALU	a	b	c
S-0470-01	S-7470-01	42mm	62mm	170mm
S-0470-02	S-7470-02	73mm	65mm	170mm



FINOCHIETTO

for adults, with suture separator

S/S	ALU	a	b	c
S-0472-01	S-7472-01	42mm	62mm	170mm





FINOCHIETTO

for adults, with 2 blades attached to each arm

S/S	ALU	a	b	c
S-0480-01	S-7480-01	10mm	15mm	120mm
S-0480-02	S-7480-02	14mm	14mm	120mm
S-0480-03	S-7480-03	14mm	20mm	120mm

THORACIC & LUNG SURGERY

Scheme of dimensions




FINOCHIETTO

3 pairs of exchangeable blades


S/S	ALU
S-0490-00	S-7490-00
blades	
S-0490-30	S-7490-30
S-0490-40	S-7490-40
S-0490-50	S-7490-50
frame only 150mm	
S-0490-01	S-7490-01

S/S	ALU
S-0491-00	S-7491-00
blades	
S-0491-40	S-7491-40
S-0491-50	S-7491-50
S-0491-60	S-7491-60
frame only 250mm	
S-0491-01	S-7491-01


THORACIC & LUNG SURGERY




S/S Alu
S-0490-30 S-7490-30
30x30mm




S/S Alu
S-0490-40 S-7490-40
40x40mm




S/S Alu
S-0490-50 S-7490-50
50x50mm



S/S Alu
S-0491-40 S-7491-40
40x100mm




S/S Alu
S-0491-50 S-7491-50
50x80mm




S/S Alu
S-0491-60 S-7491-60
60x60mm


Optional Exchangeable Blades



S/S Alu
S-0491-70 S-7491-70
70x50mm




S/S Alu
S-0491-80 S-7491-80
80x60mm



S/S Alu
S-0491-90 S-7491-90
80x90mm

THORACIC & LUNG SURGERY

Scheme of dimensions



BURFORD

for children
with 2 pairs of blades

S/S	ALU	a	b	c
S-0500-01	S-7500-01	30mm 41mm	42mm 42mm	130mm
S-0500-02	S-7500-02	34mm 48mm	42mm 42mm	135mm





BURFORD

for adults
with 2 pairs of blades

S/S	ALU	a	b	c
S-0502-01	S-7502-01	47mm 65mm	62mm 62mm	200mm
S-0502-02	S-7502-02	47mm 65mm	62mm 62mm	260mm
S-0502-03	S-7502-03	47mm 65mm	62mm 62mm	300mm


THORACIC & LUNG SURGERY



THORACIC & LUNG SURGERY




Scheme of dimensions



MORSE STERNAL



S/S	ALU	a	b	c
S-0520-01	S-7520-01	30mm	35mm	200mm



BRAASTAD


S/S	ALU
S-0530-01	S-7530-01
with 2 pairs of blades	
S-0530-65	S-7530-65
S-0530-10	S-7530-10

THORACIC & LUNG SURGERY




S/S	ALU	a	b	c
S-0540-01	36.7540.01	20mm	20mm	160mm

H-ANKENEY
child size, with double swivel blades
distance between blades 30mm, suture separator



NISSEN			
S/S	a	b	c
S-0550-01	45mm	75mm	230mm
	40mm	53mm	
S-0550-02	45mm	75mm	250mm
	40mm	53mm	





ANKENEY			
S/S	a	b	c
S-0560-01	52mm	35mm	178mm

52x35mm


S-0560-06
3 pairs of blades
52x35mm

THORACIC & LUNG SURGERY




TUFFIER

S/S	a	b	c
S-0570-14	32mm	19mm	110mm



TUFFIER


S/S	a	b	c
S-0572-17	43mm	53mm	130mm




TUFFIER

S/S	a	b	c
S-0574-16	50mm	45mm	165mm

THORACIC & LUNG SURGERY



HAIGHT			
S/S	a	b	c
S-0580-14	30mm	30mm	90mm




HAIGHT			
S/S	a	b	c
S-0582-15	30mm	30mm	115mm




HERTZLER			
S/S	a	b	c
S-0590-15	20mm	25mm	105mm


THORACIC & LUNG SURGERY



THORACIC & LUNG SURGERY



THORACIC & LUNG SURGERY



DE BAKEY

S-0640-00

with 3 pairs of blades

S-0640-01 / S-0640-02 / S-0640-03


S-0640-22 spreader only



Blades
S-0640-01
30 x 30 mm




Blades
S-0640-02
40 x 40 mm



Blades
S-0640-03
50 x 50 mm


THORACIC & LUNG SURGERY




DE BAKEY

S-0642-00
 with 3 pairs of blades
 S-0642-01/S-0642-02/S-0642-03
 S-0642-28 spreader only


blades
 S-0642-03
 40 x 100 mm




blades
 S-0642-01
 60 x 60 mm



blades
 S-0642-04
 80 x 60 mm




blades
 S-0642-02
 50 x 80 mm



blades
 S-0642-05
 50 x 66 mm



THORACIC & LUNG SURGERY




QUERVAIN
S-0650-00
35 X 50 X 180 mm

THORACIC & LUNG SURGERY




42 x 65 mm



S-0660-65

75 x 75 mm



S-0660-75

S-0660-00 complete consisting of:
S-0660-01 retractor only
S-0660-65 1 pair blades 42x65mm
S-0660-75 1 pair blades 75x75mm

THORACIC & LUNG SURGERY



ALLISON
S-0690-25 25.0cm
S-0690-30 30.0cm
S-0690-32 32.0cm



S-0700-26
26.0cm



CRAFOORD
S-0710-28
28.0cm



BAILEY-BABY
S-0720-14
14.0cm



BAILEY
S-0722-16 16.0cm
S-0722-20 20.0cm



BAILY-GIBBON
S-0724-16 16.0cm
S-0724-20 20.0cm



S-0730-19
19.0cm



LEMMON
S-0740-14 14.5cm
S-0740-18 18.0cm

THORACIC & LUNG SURGERY

